



Medical mounting arms and mounting systems



PREMIUM QUALITY WITH EXCEPTIONAL FLEXIBILITY

TECHNICAL DETAILS

rm-port	single	double	double
Lengths	250 mm / 325 mm	250 mm + 250 mm	325 mm + 325 mm
Monitor weight	up to 18 kg / 39.6 lbs	up to 18 kg / 39.6 lbs	up to 18 kg / 39.6 lbs
Load capacity	up to 23 kg / 50.6 lbs	up to 23 kg / 50.6 lbs	up to 18 kg / 39.6 lbs

mf-port	single	double
Lengths	350 mm	350 mm + 350 mm
Monitor weight	up to 18 kg / 39.6 lbs	up to 18 kg / 39.6 lbs
Load capacity	up to 23 kg / 50.6 lbs	up to 23 kg / 50.6 lbs

lf-port	support arm	support arms on pivot arm	
Lengths	225 mm	225 mm + 250 mm	225 mm + 325 mm
Monitor weight	up to 8 kg / 17.6 lbs	up to 8 kg / 17.6 lbs	up to 8 kg / 17.6 lbs
Load capacity	up to 8 kg / 17.6 lbs	up to 18 kg / 39.6 lbs	up to 18 kg / 39.6 lbs

flexion-port	support arms	support arms on pivot arm	
Lengths	300 mm / 400 mm	300 mm + 350 mm	400 mm + 350 mm
Monitor weight	up to 20 kg / 44 lbs	up to 14 kg / 30.8 lbs	up to 14 kg / 30.8 lbs
Load capacity	up to 20 kg / 44 lbs	up to 14 kg / 30.8 lbs	up to 14 kg / 30.8 lbs

Get in touch »

Europe: +49 89 614425-0 | sales@itd-cart.com
 America: +1 800 947-3901 | salesusa@itd-cart.com
 China: +86 21 5031 8818 | saleschina@itd-cart.com
 India: +91 86 0015 4064 | salesindia@itd-cart.com
 Australia: +61 3 9330 5050 | salesaustralia@itd-cart.com

www.itd-cart.com




MOUNTING SYSTEMS

Medical mounting system modul-port


Different iTD profiles – in lengths from 320 mm to 1920 mm




A wide range of monitor mounts – to accommodate standard monitors



Large portfolio of accessories and components – for various applications



Mechanical / electrical height adjustment



DS-0292.799 – 01/2025. Printed in Germany. All rights to make changes reserved.

